

## Anthocyanin Assay Kit (Microanalysis)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years  
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Anthocyanins are a class of natural pigments that are easily soluble in polar solvents and belong to flavonoids. They widely exist in the roots, stems, leaves, flowers and fruits of plants, giving them different colors from red to purple, and are the main coloring substances in plants.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481