

## Anthocyanin Assay Kit (Spectrophotometry)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years  
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Anthocyanins are a class of naturally occurring pigments that are easily soluble in polar solvents and belong to the flavonoid family. They are widely distributed in the roots, stems, leaves, flowers, and fruits of plants, giving them colors ranging from red to purple. Anthocyanins are the primary pigments responsible for the coloration of plants.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481