

MDL-72527 free base

Chemical Properties


CAS No. : 99207-33-7

Formula: C₁₂H₂₀N₂

Molecular Weight: 192.3

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	MDL-72527 free base is a polyamine oxidase inactivator.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2002 mL	26.001 mL	52.0021 mL
5 mM	1.040 mL	5.2002 mL	10.4004 mL
10 mM	0.520 mL	2.6001 mL	5.2002 mL
50 mM	0.104 mL	0.520 mL	1.040 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481