

Goat Anti-Human IgG, Fcγ fragment specific Secondary Antibody-FITC

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	FITC
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

Applications

Application:	FCM
Recommended	FCM:1:50-1:200.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Antigen Details

Antigen Species:	Human
------------------	-------

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481