

3-hydroxyandrostane-6,17-dione

Chemical Properties


CAS No. :

Formula: C₁₉H₂₈O₃

Molecular Weight: 304.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-hydroxyandrostane-6,17-dione is a useful organic compound for research related to life sciences and the catalog number is T131625.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2848 mL	16.4241 mL	32.8483 mL
5 mM	0.657 mL	3.2848 mL	6.5697 mL
10 mM	0.3285 mL	1.6424 mL	3.2848 mL
50 mM	0.0657 mL	0.3285 mL	0.657 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481