

Compound N094-0010

Chemical Properties


CAS No. :

Formula: C13H4F5NO4S

Molecular Weight: 365.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N094-0010 is a useful organic compound for research related to life sciences and the catalog number is T131161.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.738 mL	13.690 mL	27.380 mL
5 mM	0.5476 mL	2.738 mL	5.476 mL
10 mM	0.2738 mL	1.369 mL	2.738 mL
50 mM	0.0548 mL	0.2738 mL	0.5476 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481