


## Chloronectrin

## Chemical Properties

CAS No. :	38965-84-3
Formula:	C <sub>25</sub> H <sub>33</sub> ClO <sub>6</sub>
Molecular Weight:	464.979
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Chloronectrin exhibits activity against Gram-positive bacteria.
Targets(IC50)	Antibacterial,Antibiotic

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1506 mL	10.7532 mL	21.5063 mL
5 mM	0.4301 mL	2.1506 mL	4.3013 mL
10 mM	0.2151 mL	1.0753 mL	2.1506 mL
50 mM	0.043 mL	0.2151 mL	0.4301 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481