


## Feigrisolide B

## Chemical Properties

CAS No. :	313949-17-6
Formula:	C <sub>11</sub> H <sub>20</sub> O <sub>4</sub>
Molecular Weight:	216.274
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Feigrisolide B is a lactone produced by the gray Streptomyces. It exhibits antifungal activity against Trichosporon pullulans and inhibitory effects on Coxsackie virus B3. Additionally, Feigrisolide B has moderate inhibitory activity against 3 $\alpha$ -hydroxysteroid dehydrogenase.
Targets(IC50)	Antibiotic, Antifungal

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6238 mL	23.1192 mL	46.2385 mL
5 mM	0.9248 mL	4.6238 mL	9.2477 mL
10 mM	0.4624 mL	2.3119 mL	4.6238 mL
50 mM	0.0925 mL	0.4624 mL	0.9248 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481