


## o-Diaminobenzene-D4

## Chemical Properties

CAS No. :	291765-93-0
Formula:	C <sub>6</sub> H <sub>8</sub> N <sub>2</sub>
Molecular Weight:	112.166
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	o-Diaminobenzene-D4 is the deuterium-labeled version of o-Diaminobenzene.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.915 mL	44.5752 mL	89.1504 mL
5 mM	1.783 mL	8.915 mL	17.8301 mL
10 mM	0.8915 mL	4.4575 mL	8.915 mL
50 mM	0.1783 mL	0.8915 mL	1.783 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481