


## Panclcin B

## Chemical Properties

CAS No. :	160700-42-5
Formula:	C <sub>26</sub> H <sub>47</sub> N <sub>0</sub> O <sub>5</sub>
Molecular Weight:	453.655
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Panclcin B is an effective inhibitor of pancreatic lipase (pancreaticlipase).
Targets(IC50)	Lipase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2043 mL	11.0215 mL	22.0429 mL
5 mM	0.4409 mL	2.2043 mL	4.4086 mL
10 mM	0.2204 mL	1.1021 mL	2.2043 mL
50 mM	0.0441 mL	0.2204 mL	0.4409 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481