


CDK1/2/4-IN-2

## Chemical Properties

CAS No. : 677701-62-1  
 Formula: C<sub>18</sub>H<sub>16</sub>N<sub>4</sub>O  
 Molecular Weight: 304.346  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	WAY-643018 (Compound 2) is an indolecarboxamide derivative that functions as an inhibitor of CDK1, CDK2, and CDK4. CDK1/2/4-IN-2 shows potential for use in cancer research.
Targets(IC50)	CDK

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2857 mL	16.4285 mL	32.8569 mL
5 mM	0.6571 mL	3.2857 mL	6.5714 mL
10 mM	0.3286 mL	1.6428 mL	3.2857 mL
50 mM	0.0657 mL	0.3286 mL	0.6571 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481