


KYS-1

## Chemical Properties

CAS No. : 592544-81-5  
Formula: C<sub>61</sub>H<sub>82</sub>N<sub>16</sub>O<sub>18</sub>  
Molecular Weight: 1327.4  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	KYS-1 is a selective substrate for cathepsin E, with a Km value of 1.9 μM. It is incubated with cathepsin E (1.0 ng) at an optimal pH of approximately 4.0.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7534 mL	3.7668 mL	7.5335 mL
5 mM	0.1507 mL	0.7534 mL	1.5067 mL
10 mM	0.0753 mL	0.3767 mL	0.7534 mL
50 mM	0.0151 mL	0.0753 mL	0.1507 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481