


Siamine

Chemical Properties

CAS No. :	60352-12-7
Formula:	C ₁₀ H ₉ NO ₃
Molecular Weight:	191.183
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Siamine is an isoquinolone alkaloid that can be extracted from Cassia siamea and serves as an antioxidant.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2307 mL	26.1534 mL	52.3067 mL
5 mM	1.0461 mL	5.2307 mL	10.4613 mL
10 mM	0.5231 mL	2.6153 mL	5.2307 mL
50 mM	0.1046 mL	0.5231 mL	1.0461 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481