


Napyradiomycin A2

Chemical Properties

CAS No. :	111216-62-7
Formula:	C ₂₅ H ₃₀ Cl ₂ O ₆
Molecular Weight:	497.408
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Napyradiomycin A2 is an antibiotic that exhibits activity against Gram-positive bacteria and mycobacteria.
Targets(IC50)	Antibacterial, Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0104 mL	10.0521 mL	20.1041 mL
5 mM	0.4021 mL	2.0104 mL	4.0208 mL
10 mM	0.201 mL	1.0052 mL	2.0104 mL
50 mM	0.0402 mL	0.201 mL	0.4021 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481