


Galacardin A

Chemical Properties

CAS No. :	137801-55-9
Formula:	C ₁₀₁ H ₁₂₁ Cl ₂ N ₉ O ₄₆
Molecular Weight:	2267.98
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Galacardin A exhibits activity against Gram-positive bacteria.
Targets(IC50)	Antibacterial, Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4409 mL	2.2046 mL	4.4092 mL
5 mM	0.0882 mL	0.4409 mL	0.8818 mL
10 mM	0.0441 mL	0.2205 mL	0.4409 mL
50 mM	0.0088 mL	0.0441 mL	0.0882 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481