


Elloramycin

Chemical Properties

CAS No. :	97218-42-3
Formula:	C32H36O15
Molecular Weight:	660.619
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Elloramycin exhibits mild activity against Gram-positive bacteria, streptomycetes, and L-1210 leukemia cells, but it is ineffective against P388 leukemia cells in mice. Additionally, it demonstrates significant inhibitory effects on streptomycetes, including actinomyces, within the brackets.
Targets(IC50)	Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5137 mL	7.5686 mL	15.1373 mL
5 mM	0.3027 mL	1.5137 mL	3.0275 mL
10 mM	0.1514 mL	0.7569 mL	1.5137 mL
50 mM	0.0303 mL	0.1514 mL	0.3027 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481