


## Collinone

## Chemical Properties

CAS No. :	343315-22-0
Formula:	C <sub>27</sub> H <sub>18</sub> O <sub>12</sub>
Molecular Weight:	534.425
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Collinone is a recombinant polyketide antibiotic with antimicrobial activity against Gram-positive bacteria and exhibits cytotoxicity.
Targets(IC50)	Antibacterial, Antibiotic

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8712 mL	9.3558 mL	18.7115 mL
5 mM	0.3742 mL	1.8712 mL	3.7423 mL
10 mM	0.1871 mL	0.9356 mL	1.8712 mL
50 mM	0.0374 mL	0.1871 mL	0.3742 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481