


(S)-Albuterol

Chemical Properties

CAS No. :	34271-50-6
Formula:	C ₁₃ H ₂₁ NO ₃
Molecular Weight:	239.311
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-Albuterol functions as a muscarinic receptor and phospholipase C activator, enhancing the intracellular free calcium concentration in airway smooth muscle cells.
Targets(IC50)	Calcium Channel,AChR,Phospholipase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1787 mL	20.8934 mL	41.7868 mL
5 mM	0.8357 mL	4.1787 mL	8.3574 mL
10 mM	0.4179 mL	2.0893 mL	4.1787 mL
50 mM	0.0836 mL	0.4179 mL	0.8357 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481