


Antiviral agent 67

Chemical Properties

CAS No. : 25111-96-0
 Formula: C₁₉H₁₉N₃O
 Molecular Weight: 305.374
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Antiviralagent 67 (compound PC6) is an inhibitor of DENVNS5 (RNA-dependent RNA polymerase) with a K_i value of 1.12 nM.
Targets(IC50)	DNA/RNA Synthesis

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2747 mL	16.3736 mL	32.7472 mL
5 mM	0.6549 mL	3.2747 mL	6.5494 mL
10 mM	0.3275 mL	1.6374 mL	3.2747 mL
50 mM	0.0655 mL	0.3275 mL	0.6549 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481