


Hexamidine

Chemical Properties

CAS No. :	3811-75-4
Formula:	C ₂₀ H ₂₆ N ₄ O ₂
Molecular Weight:	354.446
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hexamidine exhibits amoebicidal activity and is applicable in the study of Acanthamoeba keratitis (Acanthamoebakeratitis).
Targets(IC50)	Parasite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8213 mL	14.1064 mL	28.2127 mL
5 mM	0.5643 mL	2.8213 mL	5.6425 mL
10 mM	0.2821 mL	1.4106 mL	2.8213 mL
50 mM	0.0564 mL	0.2821 mL	0.5643 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481