


Secutrelvir

Chemical Properties

CAS No. :	2996148-73-1
Formula:	C ₂₃ H ₁₆ Cl ₂ F ₃ N ₅ O ₂
Molecular Weight:	522.307
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Secutrelvir inhibits the 3CL protease (3CL protease), thereby exhibiting antiviral activity by preventing the replication of SARS-CoV-2.
Targets(IC50)	SARS-CoV, Virus Protease

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9146 mL	9.5729 mL	19.1457 mL
5 mM	0.3829 mL	1.9146 mL	3.8291 mL
10 mM	0.1915 mL	0.9573 mL	1.9146 mL
50 mM	0.0383 mL	0.1915 mL	0.3829 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481