


10-Bromodecanoic acid

Chemical Properties

CAS No. :	50530-12-6
Formula:	C ₁₀ H ₁₉ BrO ₂
Molecular Weight:	251.161
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	10-Bromodecanoic acid (Decanoic acid, 10-bromo-) is a linker used in the synthesis of PROTACs, such as PROTAC KDM4 degrader-1.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9815 mL	19.9076 mL	39.8153 mL
5 mM	0.7963 mL	3.9815 mL	7.9631 mL
10 mM	0.3982 mL	1.9908 mL	3.9815 mL
50 mM	0.0796 mL	0.3982 mL	0.7963 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481