

E3 Ligase Ligand-linker Conjugate 176

Chemical Properties

CAS No. : 3005273-99-1


Formula: C₁₉H₂₅N₅O₂

Molecular Weight: 355.434

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	E3LigaseLigand-linker Conjugate 176 is an E3 ubiquitin ligase ligand-linker conjugate (E3ligaseligand-linker conjugate) that can be utilized in the synthesis of BMS-986458.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8135 mL	14.0675 mL	28.1349 mL
5 mM	0.5627 mL	2.8135 mL	5.627 mL
10 mM	0.2813 mL	1.4067 mL	2.8135 mL
50 mM	0.0563 mL	0.2813 mL	0.5627 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481