

Mal-amido-PEG8-val-gly

Chemical Properties


CAS No. : 2353409-54-6

Formula: C₃₃H₅₆N₄O₁₅

Molecular Weight: 748.815

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Mal-amido-PEG8-Val-Gly is an ADC linker containing a maleimide moiety.
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3354 mL	6.6772 mL	13.3543 mL
5 mM	0.2671 mL	1.3354 mL	2.6709 mL
10 mM	0.1335 mL	0.6677 mL	1.3354 mL
50 mM	0.0267 mL	0.1335 mL	0.2671 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481