


1-Linoleoyl-(2S)-glycerol

Chemical Properties

CAS No. : 67968-46-1
 Formula: C₂₁H₃₈O₄
 Molecular Weight: 354.524
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Linoleoyl-(2S)-glycerol is a lipid structurally similar to linoleic acid. The hydroxyl group can undergo further derivatization or be replaced by other active functional
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8207 mL	14.1036 mL	28.2072 mL
5 mM	0.5641 mL	2.8207 mL	5.6414 mL
10 mM	0.2821 mL	1.4104 mL	2.8207 mL
50 mM	0.0564 mL	0.2821 mL	0.5641 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481