


SU-740

Chemical Properties

CAS No. : 134336-72-4
 Formula: C₂₆H₃₀O₅
 Molecular Weight: 422.513
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	SU-740 is a chalcone derivative with anti-ulcer activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3668 mL	11.834 mL	23.6681 mL
5 mM	0.4734 mL	2.3668 mL	4.7336 mL
10 mM	0.2367 mL	1.1834 mL	2.3668 mL
50 mM	0.0473 mL	0.2367 mL	0.4734 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481