


## FAM azide, 5-isomer

## Chemical Properties

CAS No. :	510758-23-3
Formula:	C <sub>24</sub> H <sub>18</sub> N <sub>4</sub> O <sub>6</sub>
Molecular Weight:	458.423
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	FAM azide, 5-isomer, is a reagent used for Click Chemistry labeling.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1814 mL	10.907 mL	21.8141 mL
5 mM	0.4363 mL	2.1814 mL	4.3628 mL
10 mM	0.2181 mL	1.0907 mL	2.1814 mL
50 mM	0.0436 mL	0.2181 mL	0.4363 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481