

N-Boc-Tris

Chemical Properties

CAS No. :	146651-71-0
Formula:	C ₉ H ₁₉ NO ₅
Molecular Weight:	221.251
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	N-Boc-Tris is a protective reagent used in organic synthesis and pharmaceutical development to shield amines with the Boc (tert-butoxycarbonyl) group, thus preventing unnecessary side reactions.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5198 mL	22.5989 mL	45.1977 mL
5 mM	0.904 mL	4.5198 mL	9.0395 mL
10 mM	0.452 mL	2.2599 mL	4.5198 mL
50 mM	0.0904 mL	0.452 mL	0.904 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481