


2,2'-(Piperazine-1,4-diyl)diethanamine

Chemical Properties

CAS No. :	6531-38-0
Formula:	C ₈ H ₂₀ N ₄
Molecular Weight:	172.271
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,2'-(Piperazine-1,4-diyl)diethanamine is a piperazine compound characterized by the presence of two N-ethylamine groups. Such molecules can be utilized in pharmaceutical synthesis or in the design of innovative ionizable amino lipid head groups for lipid nanoparticles.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8048 mL	29.0242 mL	58.0484 mL
5 mM	1.161 mL	5.8048 mL	11.6097 mL
10 mM	0.5805 mL	2.9024 mL	5.8048 mL
50 mM	0.1161 mL	0.5805 mL	1.161 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481