


TTBK1/2-IN-2

## Chemical Properties

CAS No. : 2857982-35-3  
 Formula: C<sub>22</sub>H<sub>24</sub>N<sub>4</sub>O  
 Molecular Weight: 360.452  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	TTBK1/2-IN-2 (compound 9) is a potent inhibitor of Tau tubulin kinase 1 (TTBK1) and TTBK2, with IC <sub>50</sub> values of 384 nM and 175 nM, respectively. In human induced pluripotent stem cells (iPSC), TTBK1/2-IN-2 demonstrates a significant ciliogenesis phenotype.
Targets(IC <sub>50</sub> )	Microtubule Associated

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7743 mL	13.8715 mL	27.7431 mL
5 mM	0.5549 mL	2.7743 mL	5.5486 mL
10 mM	0.2774 mL	1.3872 mL	2.7743 mL
50 mM	0.0555 mL	0.2774 mL	0.5549 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481