


MAGL Ligand-Linker Conjugate 1

Chemical Properties

CAS No. :	2725023-11-8
Formula:	C ₂₈ H ₂₈ N ₂ O ₉
Molecular Weight:	536.53
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	MAGL Ligand-Linker Conjugate 1 is a compound that combines a MAGL ligand with a linker. This conjugate serves as an essential component for synthesizing the PROTAC degrader, known as PROTACMAGLdegrader-1.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8638 mL	9.3191 mL	18.6383 mL
5 mM	0.3728 mL	1.8638 mL	3.7277 mL
10 mM	0.1864 mL	0.9319 mL	1.8638 mL
50 mM	0.0373 mL	0.1864 mL	0.3728 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481