


## Bagremycin A

## Chemical Properties

CAS No. :	377775-45-6
Formula:	C <sub>15</sub> H <sub>13</sub> NO <sub>3</sub>
Molecular Weight:	255.269
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Bagremycin A is identified in the actinobacterium strain Streptomyces (Streptomyces) Tu 4128. It exhibits weak antimicrobial activity against Gram-positive bacteria, Saccharomyces cerevisiae, and Candida albicans.
Targets(IC50)	Antibacterial, Antifungal

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9174 mL	19.5871 mL	39.1742 mL
5 mM	0.7835 mL	3.9174 mL	7.8348 mL
10 mM	0.3917 mL	1.9587 mL	3.9174 mL
50 mM	0.0783 mL	0.3917 mL	0.7835 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481