


PCSK9 ligand 2

Chemical Properties

CAS No. :	3039841-82-9
Formula:	C16H19N3O3S
Molecular Weight:	333.405
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	PCSK9 ligand 2 is the target protein ligand of PROTAC PCSK9 autophagic degrader 2, and it is applicable for research related to anti-atherosclerosis.
Targets(IC50)	Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9993 mL	14.9966 mL	29.9931 mL
5 mM	0.5999 mL	2.9993 mL	5.9986 mL
10 mM	0.2999 mL	1.4997 mL	2.9993 mL
50 mM	0.060 mL	0.2999 mL	0.5999 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481