


(-)-Codonopsine

Chemical Properties

CAS No. :	26989-20-8
Formula:	C ₁₄ H ₂₁ NO ₄
Molecular Weight:	267.321
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(-)-Codonopsine is an alkaloid found in <i>Codonopsis clematidea</i> , exhibiting antibacterial and hypotensive properties. The minimum inhibitory concentration (MIC) value against <i>Escherichia coli</i> (<i>E. coli</i>) is 64 µg/mL. (-)-Codonopsine shows promise for research in anti-infective and cardiovascular applications.
Targets(IC50)	Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7408 mL	18.7042 mL	37.4083 mL
5 mM	0.7482 mL	3.7408 mL	7.4817 mL
10 mM	0.3741 mL	1.8704 mL	3.7408 mL
50 mM	0.0748 mL	0.3741 mL	0.7482 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481