

E3 Ligase Ligand-linker Conjugate 183

Chemical Properties

CAS No. : 3081572-35-9


Formula: C₂₃H₂₃Br₃INO₅

Molecular Weight: 760.049

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	E3LigaseLigand-linker Conjugate 183 is an E3 ligase ligand-linker conjugate of PCSK9 autophagy degrader 2. It comprises the LC3B ligand 2 and its corresponding linker.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3157 mL	6.5785 mL	13.157 mL
5 mM	0.2631 mL	1.3157 mL	2.6314 mL
10 mM	0.1316 mL	0.6579 mL	1.3157 mL
50 mM	0.0263 mL	0.1316 mL	0.2631 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481