


## Thionordiazepam

## Chemical Properties

CAS No. : 4547-02-8  
Formula: C<sub>15</sub>H<sub>11</sub>ClN<sub>2</sub>S  
Molecular Weight: 286.779  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Thionordiazepam is classified as a benzodiazepine compound.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.487 mL	17.435 mL	34.8699 mL
5 mM	0.6974 mL	3.487 mL	6.974 mL
10 mM	0.3487 mL	1.7435 mL	3.487 mL
50 mM	0.0697 mL	0.3487 mL	0.6974 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481