


Xanthinol

Chemical Properties

CAS No. : 2530-97-4
 Formula: C₁₃H₂₁N₅O₄
 Molecular Weight: 311.337
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Xanthinol is an effective modulator of tumor hemodynamics.
Targets(IC50)	PDGFR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2119 mL	16.0596 mL	32.1192 mL
5 mM	0.6424 mL	3.2119 mL	6.4238 mL
10 mM	0.3212 mL	1.606 mL	3.2119 mL
50 mM	0.0642 mL	0.3212 mL	0.6424 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481