


FAP targeting peptide-PEG2 conjugate

Chemical Properties

CAS No. :	3076328-71-4
Formula:	C73H102N16O18
Molecular Weight:	1491.69
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	FAP-targeting peptide-PEG2 conjugate (Example A1) is the polypeptide and linker conjugate component of Unlabeled FXX489. Unlabeled FXX489 serves as a ligand targeting fibroblast activation protein (FAP).
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6704 mL	3.3519 mL	6.7038 mL
5 mM	0.1341 mL	0.6704 mL	1.3408 mL
10 mM	0.067 mL	0.3352 mL	0.6704 mL
50 mM	0.0134 mL	0.067 mL	0.1341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481