


## Isopropyl benzenesulfonate

## Chemical Properties

CAS No. : 6214-18-2  
Formula: C<sub>9</sub>H<sub>12</sub>O<sub>3</sub>S  
Molecular Weight: 200.255  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Isopropyl benzenesulfonate (Isopropyl besylate) is a sulfonate ester that is also considered a potential genotoxic impurity.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9935 mL	24.9675 mL	49.9351 mL
5 mM	0.9987 mL	4.9935 mL	9.987 mL
10 mM	0.4994 mL	2.4968 mL	4.9935 mL
50 mM	0.0999 mL	0.4994 mL	0.9987 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481