

Tos-PEG22-Tos

Chemical Properties

CAS No. : 35164-96-6

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Tos-PEG22-Tos is a PEG linker that incorporates two tosyl (p-toluenesulfonyl) groups. The hydrophilic PEG spacer enhances solubility in aqueous media, and the tosyl group serves as an excellent leaving group in nucleophilic substitution reactions.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481