

## Keratin sulfate

### Chemical Properties

CAS No. : 169799-44-4

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Keratin sulfate is a heparin-related proteoglycan that can interact with receptors or alter their stability and function. Unlike some other compounds, keratin sulfate does not enhance the mitogenic effect of FGF-2.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481