

Poly(D,L-lactide-co-glycolide)

Chemical Properties

CAS No. : 26780-50-7

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Poly(D,L-lactide-co-glycolide) is a synthetic polymer material known for its excellent biocompatibility. It is utilized in research concerning drug delivery systems.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481