

YY2201

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	YY2201 is a potent, selective inhibitor of ATR with an IC50 of 8 nM. It exhibits over 200 times greater selectivity for ATR compared to mTOR. Additionally, YY2201 is capable of inhibiting tumor progression across a wide range of cancer types.
Targets(IC50)	ATM/ATR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481