

## Prostaglandin D synthase

## Chemical Properties

CAS No. : 65802-85-9

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Prostaglandin D synthase is a biomarker for meningioma cells and coronary artery disease. Lipocalin-type prostaglandin D synthase (L-PGDS) is found in human coronary atherosclerotic plaques and can be detected in human serum.
Targets(IC50)	PGE Synthase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481