

F11R peptide TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

F11R peptide TFA is a molecule of the F11 receptor. It inhibits platelet aggregation induced by anti-F11R antibodies and prevents platelet adhesion to endothelial cells inflamed by cytokines (TNF α and INF- γ). F11R peptide TFA is applicable in research on atherosclerotic plaques and related thrombotic diseases.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481