

PLK1-IN-14

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PLK1-IN-14 (Compound Hit-4) is a selective PLK1 inhibitor with an IC50 of 22.61 pM for PLK1, and an IC50 of 0.09 nM for inhibiting the proliferation of DU-145 prostate cancer cells. It is applicable in prostate cancer research.
Targets(IC50)	PLK

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481