

SrtA-IN-1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | SrtA-IN-1 (Compound T10) is an inhibitor of Streptococcus pyogenes sortase A (SpSrtA) with an IC50 value of 0.7 μ M. It is applicable for research in the field of anti-infective agents. |
| Targets(IC50) | Antibacterial |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481