

FITC-PEG-NHS (MW 400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FITC-PEG-NHS (MW 400) is a fluorescent dye composed of polyethylene glycol (PEG), fluorescein isothiocyanate (FITC), and N-hydroxysuccinimide (NHS) groups. The NHS group is capable of reacting with amine-containing groups. FITC-PEG-NHS (MW 400) is applicable for fluorescence labeling and imaging.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481