

DSPE-PEG1000-TH

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG1000-TH is a PEG compound composed of DSPE and a pH-responsive cell-penetrating peptide (TH). In acidic environments, such as the tumor microenvironment, TH is activated, enabling the selective delivery of small molecules, oligonucleotides, proteins, and other agents into tumor cells.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481