

EGFR-IN-155

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	EGFR-IN-155 (compound 13a) is an EGFR inhibitor with IC50 values of 0.14 nM for EGFR TK and 0.18 nM for EGFR L858R. EGFR-IN-155 suppresses tumor growth, causes cell cycle arrest at the S phase, and induces apoptosis (cell death).
Targets(IC50)	Apoptosis,EGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481